# bluephase® C8

Licence to cure

reddot design award  
honourable mention 2008


For dental use only!  
Caution: U.S. Federal law restricts this device to  
sale by or on the order of a licensed dentist  
Made in Austria  
Ivoclar Vivadent AG, FL-9494 Schaan /  
Liechtenstein

ivoclar  
vivadent®  
clinical

# bluephase® C8

LED for every use

## Technical data

Wavelength range:	380-515 nm
Light intensity:	800 mW/cm <sup>2</sup> ± 10%
Curing modes:	High Power, Soft Start, Low Power
Exposure time:	20 sec 15 sec (Tetric EvoCeram / IPS Empress Direct)
Light probe:	10 mm
Power supply:	Mains operated
Operating voltage:	100-240 VAC / 50-60 Hz
Warranty	3 years

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

## bluephase® C8

Light intensity	800 mW/cm <sup>2</sup> ± 10%
Exposure time for composites	20 sec.
Exposure time for Tetric EvoCeram / IPS Empress Direct	15 sec.

### Filling materials

#### Composite

Tetric EvoCeram / Tetric EvoFlow / IPS Empress Direct / Artemis / 4 Seasons / InTen-S / Tetric new / Tetric N-Ceram / Tetric N-Flow	15 sec.
Heliomolar / Heliomolar HB / Heliomolar Flow	20 sec.
Others <sup>1)</sup>	

#### Compomer<sup>2)</sup>

#### Compoglass F / Compoglass Flow

40 sec.

### Indirect restorations / Luting materials

#### Variolink II Base<sup>3)</sup> / Variolink Veneer

per mm ceramic: 10 sec. per segment

#### Variolink Ultra<sup>4)</sup> / Dual Cement<sup>4)</sup>

per mm ceramic: 10 sec. per segment

#### Variolink II<sup>4)</sup>

### Miscellaneous

#### Helioseal / Helioseal F / Helioseal Clear

20 sec.

#### Monopaque

30 sec.

#### MultiCore Flow / Multicore HB

40 sec.

#### CoroSeal

20 sec.

#### Heliosit Orthodontic

20 sec.

### Adhesive

#### AdheSE / AdheSE One F / Excite F / Excite F DSC / Tetric N-Bond / Tetric N-Bond Self-Etch / Heliobond (Syntac)

10 sec.

### Temporary materials

#### Systemp.link

10 sec.

#### Systemp.inlay/onlay

10 sec.

#### Fermit / Fermit N

10 sec.

### Miscellaneous

#### Heliosit Orthodontic

20 sec.

#### CoroSeal Bonding

10 sec.

#### Vivaglass Liner

20 sec.

Dear Customer,

Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The polymerization light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased bluephase C8. bluephase® C8 is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.

These Operating Instructions will help you safely to start up the device, make full use of its capabilities and ensure a long service life.

Should you have any further questions, please do not hesitate to contact us (see addresses on the back cover).

Your Ivoclar Vivadent Team

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case with dentin shades)

2) Applies to a maximum layer thickness of 3 mm

3) Applies to light-curing (exclusive use of base paste)

4) Applies to dual-curing

# Table of Contents

<b>Safety</b>	
Intended use	8
Indication	8
Safety notes	8
<b>Start-up</b>	
Handpiece holder	10
Handpiece	10
<b>Operation</b>	
Selecting the curing program and the curing time	12
Cure Memory function	12
<b>Start</b>	
Acoustic signals	14
Light intensity	14
Measuring the light intensity	14
<b>Maintenance and cleaning</b>	
Housing	15
Light probes	15
Disposal	15
<b>What if ...</b>	16
<b>Warranty / Procedure in case of repair</b>	16
<b>Product specifications</b>	
Delivery form	17
Accessories	17
Technical data	17
Operating conditions	17
Transportation and storage conditions	17
<b>Declaration of conformity</b>	56

# Safety

## Intended use

bluephase C8 is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended use also includes the observation of the notes and regulations in these Operating Instructions.

## Indication

With its "polywave" broadband spectrum, bluephase C8 is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 380-515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporary materials as well as luting materials for brackets and dental-lab materials such as those used in the fabrication of ceramic inlays.

## Contraindications

- Materials, the polymerization of which is activated outside the wavelength range of 380-515 nm (no materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.
- This unit should not be used near flammable anaesthetics or mixtures of flammable anaesthetics with air, oxygen or nitrous oxide.
- Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.
- Patients with heart or brain pacemakers can be treated with bluephase C8. Investigations have shown that there are no (electromagnetic) incompatibilities between the curing light and pacemakers.

## Safety notes

The bluephase C8 is a medical electronic device which is subject to IEC 601-1 (EN 60601-1) and EMC standards (see appendix). It complies with the relevant EU regulations.

The apparatus has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Operating Instructions have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

NRTL (Nationally Recognized Testing Laboratory)

## Usage and liability

- bluephase C8 must only be employed for the intended use. Any other uses are contraindicated. Liability cannot be accepted for damage resulting from misuse or failure to observe the Operating Instructions.
- The user is responsible for testing bluephase C8 for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of bluephase at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent (see Accessories). The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.

## Operating voltage

Before switching on, make sure that

- the voltage indicated on the rating plate complies with the local power supply and
- the unit has acquired the ambient temperature.

## **Assumption of impaired safety**

If it has to be assumed that safe operation is no longer possible, the power must be disconnected to avoid accidental operation. This may be the case, for example, if the apparatus is visibly damaged or no longer works correctly. The apparatus is only completely disconnected from the power supply if the power cord has been pulled out of the plug-in socket.

## **Eye protection**

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are generally sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device. Therefore, they should wear protective goggles (orange) that absorb light below 515 nm.

## **Heat development**

### **(burn hazard)**

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible or reversible damage. Therefore, the recommended curing times have to be observed.

Uninterrupted curing times of more than 40 seconds on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. Polymerize indirect restorations at intermittent intervals of 20 seconds each or use external cooling with an air stream.

## Start-up

### Delivery form

Check the delivery for completeness and any possible transportation damage (see Delivery form). If components are missing or damaged, immediately contact your dealer or your Service Center.

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply. The rating plate is attached to the underside of the power pack.

Connect the power cord with the power supply and the power pack.

### **Handpiece**

Before you insert the light probe, remove the protective foil on the corresponding opening of the handpiece.

Slightly turn the light probe while attaching it to the handpiece.

After that, mount the anti-glare cone on the light probe.

# Operation

Disinfect or autoclave contaminated surfaces of the curing light as well as of the light probes and anti-glare cones before each use. Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For this purpose, check the light probe for contamination and damage as well as the light intensity at regular intervals (e.g. with the Ivoclar Vivadent radiometer "bluephase meter").

## Selecting the curing program and the curing time

The curing program and the curing time can be individually set. bluephase C8 is equipped with the following 3 curing programs for the different indications. The desired curing program is selected with the program selection buttons. The display changes accordingly (see Indicators on the hand-piece). The device comes equipped with the following pre-set programs:

## Factory settings

 <b>HIGH</b> (High Power)	<b>15 seconds</b>
 <b>LOW</b> (Low Power)	<b>10 seconds</b>
 <b>SOFT</b> (Soft start)	<b>20 seconds</b>

The intended curing time is selected using the time selection buttons. Users may choose between 5, 10, 15, 20 and 30 seconds.

Observe the Instructions for Use of the material applied when selecting the curing time.

The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 9 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

# Curing program

## bluephase® C8

<b>Light intensity</b>	800 mW/cm <sup>2</sup> ± 10%
Exposure time for composites	20 sec.
Exposure time for Tetric EvoCeram / IPS Empress Direct	15 sec.
Others <sup>1)</sup>	

### HIGH POWER Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

### Filling materials

#### Composite

Tetric EvoCeram / Tetric EvoFlow / IPS Empress Direct / Arteneis / 4 Seasons / InTen-S / Tetric new / Tetric N-Ceram / Tetric N-Flow	15 sec.
Helimolar / Helimolar HB / Helimolar Flow	20 sec.
Others <sup>1)</sup>	

#### Compomer<sup>2)</sup>

Compoglass F / Compoglass Flow	40 sec.
--------------------------------	---------

#### Indirect restorations / Luting materials


VarioLink II Base <sup>3)</sup> / VarioLink Veneer	per mm ceramic: 10 sec. per segment
VarioLink Ultra <sup>4)</sup> / Dual Cement <sup>4)</sup>	per mm ceramic: 10 sec. per segment
VarioLink II <sup>4)</sup>	

#### Miscellaneous

Helioseal / Helioseal F / Helioseal Clear	20 sec.
Monopaque	30 sec.
MultiCore Flow / Multicore HB	40 sec.
CoroSeal	20 sec.
Heliosit Orthodontic	20 sec.

### LOW POWER Program

Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities

### Adhesive

AdheSE / AdheSE One F / Excite F DSC / Tetric N-Bond / Tetric N-Bond Self-Etch / Heliobond (Syntac)	10 sec.
---	---------

#### Temporary materials


Systemp.link	10 sec.
Systemp.inlay/onlay	10 sec.
Fermi / Fermi N	10 sec.

#### Miscellaneous

Heliosit Orthodontic	20 sec.
CoroSeal Bonding	10 sec.
Vivaglass Liner	20 sec.

### SOFT START Program

Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials

### Filling materials

#### Composite

Tetric EvoCeram / Tetric EvoFlow / IPS Empress Direct / Arteneis / 4 Seasons / InTen-S / Tetric new / Tetric N-Ceram / Tetric N-Flow	20 sec.
Helimolar / Helimolar HB / Helimolar Flow	25 sec.
Others <sup>1)</sup>	

#### Compomer<sup>2)</sup>

Compoglass F / Compoglass Flow	40 sec.
--------------------------------	---------

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case with dentin shades)

2) Applies to a maximum layer thickness of 3 mm      3) Applies to light-curing (exclusive use of base paste)      4) Applies to dual-curing

# Start up

## Cure Memory function

The last settings used, together with the combination of curing program and curing time, are automatically saved.

## Start

The light is switched on by means of the start button. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again. The fan is activated simultaneously to the light. Once the curing time has elapsed, the fan continues to run for a certain time to cool the apparatus.

## Acoustic signals

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Program change
- Curing time change
- Error message

If desired, the volume of the acoustic signals can be adjusted.

For this purpose, press the blue volume button (see Product overview) while the light is off and reduce the volume or turn the acoustic signals off with the left time selection button. If you want to turn the acoustic signals back on or increase the volume, press the blue volume button and then the right time selection button.

## Light intensity

The light intensity is maintained at a consistent level during operation. If the supplied 10 mm light probe is used, the light intensity has been calibrated to  $800 \text{ mW/cm}^2 \pm 10\%$ .

The use of light probe other than the one provided has a direct influence on the light intensity emitted.

In parallel-walled light probes (10 mm), the diameter is equal at both ends. In focussing light probes ( $10 > 8 \text{ mm}$  light probe, Pin-Point light probe  $6 > 2 \text{ mm}$ ), the diameter of the rear end is larger than that of the light emission window. The incident blue light is thus bundled to a smaller surface, which increases the light intensity emitted.

Pin-Point light probes are suitable for the polymerization of confined areas, such as the attachment of veneers prior to excess removal. For thorough curing, it is necessary to change the light probe.

## Measuring the light intensity

The light intensity of the bluephase C8 and the enclosed 10 mm light probe can be checked by means of the bluephase meter.

If the measured value does not correspond with the expected light intensity, proceed as follows:

- Check the selected curing program
- Clean light sensor if it is contaminated
- Remove the light probe and clean the light emission window of the handpiece with a cotton swab dipped in alcohol.
- Clean light probe if it is contaminated (see Maintenance and cleaning)
- Replace a damaged light probe with a new one

If these measures do not improve the results, please contact your dealer or your local Service Center.

## Maintenance and cleaning

For reasons of hygiene, the use of a disposable protective sleeve for each patient is recommended. Contaminated surfaces of the device as well as light probes and anti-glare cones must be disinfected (e.g. FD333, FD366/Dürr Dental, Incidin Liquid/Ecolab) or autoclaved before each use. Make sure that no liquids or other foreign substances enter the handpiece, the handpiece holder and particularly the power pack during cleaning (risk of electrical shock). Disconnect the handpiece from the power source when cleaning it.

### Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soap solution.

### Light probes

Clean the light probe with a disinfectant and a soft cloth. Carefully remove any possible material residue, such as composite material, from the light probe with a fingernail or a plastic spatula, for instance. Do not use sharp or pointed instruments, since they may scratch the surface of the light probe and thus reduce light transmission. Light probes and anti-glare cones may be autoclaved (e.g. 134 °C, 3 bar, 5 minutes). Check light probes for damage. If a light probe is held against a light source and certain segments appear black, then the glass fibres are broken. In this case, the light probe has to be replaced with a new one.

### Disposal

The apparatus must not be disposed of in normal domestic waste. Please visit your local Ivoclar Vivadent website for detailed information about the disposal of the unit. Dispose of unserviceable polymerization lights according to the corresponding legal requirements in your country.

## What if ...

Symbol	Causes	Error rectification
	Apparatus is overheated or under-cooled	Allow the apparatus to cool down (or assume room temperature if it is undercooled) and try again after a certain time. If the error persists, please contact your dealer or your local Service Center.
	Electronic component of the power pack defective	Disconnect power pack from the power supply and then connect it again. If the error reappears, please contact your dealer or your local Service Center.
	Electronic component of the hand-piece defective	Please contact your dealer or your local Service Center.

## Warranty / Procedure in case of repair

The warranty for bluephase C8 is three years from the date of purchase. The warranty includes an LED operating life of 300 hours. Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- damage resulting from improper handling.
- damage of components resulting from wear under standard operating conditions.
- damage resulting from external influences, e.g. blows, drop to the floor.
- damage resulting from incorrect set-up or installation.

- damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.
- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, handpiece holder, power cord and light probe) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Center. In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on the back cover). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

## Product specifications

### Delivery form

1 Handpiece holder,  
1 Handpiece including power cord and power pack,  
1 10 mm light probe,  
3 Anti-glare cones,  
1 Set of Operating Instructions

### Accessories

The following accessories are available for bluephase C8:

REF.	Description
608537	Light probe 10 mm, black (G2)
608538	Pin-point light probe 6/2 mm, black (G2)
551756	Anti-glare cones
592496	Anti-glare shield
627436	Handpiece holder C8 (G2)
607922	bluephase meter sleeves
608554	

### Technical data

Operating voltage	5 VDC
Power Pack	100–240 VAC / 50-60 Hz / max. 0.4 A Output 5 VDC TYP 15.2630 Manufacturer: Friwo

### Operating conditions

Temperature	+10 °C to +40 °C
Relative humidity	30 % to 75 %
Ambient pressure	700 hPa to 1060 hPa
Dimensions of handpiece holder	L=205mm; B=150mm; H=85mm
Weight of the handpiece holder	190 g
Light source	polywave LED
Wavelength range	380–515 nm
Light intensity	800 mW/cm <sup>2</sup> ± 10 %
Operating	6 min on / 4 min off (intermittently)
Dimensions of Handpiece:	L=260mm; B=42mm; H=120mm
Weight of the handpiece	200 g

### Transportation and storage conditions:

Temperature	-20 °C to +60 °C (-4 to 140 °F)
Relative humidity	10 % to 75 %
Ambient pressure	500 hPa to 1060 hPa

bluephase C8 has to be stored in closed, roofed rooms.  
Protect the device from severe jarring.

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# Appendix

## Guidance and manufacturer's declaration - electromagnetic emission

The bluephase C8 is intended for use in the electromagnetic environment specified below. The customer or the user of the bluephase C8 should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The bluephase C8 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The bluephase C8 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	P < 75W (pass without test)
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	P < 75W (pass without test)

## Guidance and manufacturer's declaration - electromagnetic immunity

The bluephase C8 is intended for use in the electromagnetic environment specified below. The customer or the user of the bluephase C8 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U <sub>r</sub> (>95 % dip in U <sub>r</sub> ) for 0.5 cycle  40 % U <sub>r</sub> (60 % dip in U <sub>r</sub> ) for 5 cycles  70 % U <sub>r</sub> (30 % dip in U <sub>r</sub> ) for 25 cycles  <5 % U <sub>r</sub> (>95 % dip in U <sub>r</sub> ) for 5 sec	<5 % U <sub>r</sub> (>95 % dip in U <sub>r</sub> ) for 0.5 cycle  40 % U <sub>r</sub> (60 % dip in U <sub>r</sub> ) for 5 cycles  70 % U <sub>r</sub> (30 % dip in U <sub>r</sub> ) for 25 cycles  <5 % U <sub>r</sub> (>95 % dip in U <sub>r</sub> ) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of the bluephase C8 requires continued operation during power mains interruptions, it is recommended that the bluephase C8 be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE: U<sub>r</sub> is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic immunity

The bluephase C8 is intended for use in the electromagnetic environment specified below. The customer or the user of the bluephase C8 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the bluephase C8, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			<b>Recommended separation distance</b>
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the bluephase C8 is used exceeds the applicable RF compliance level above, the bluephase C8 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the bluephase C8.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the bluephase C8

The bluephase C8 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the bluephase C8 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the bluephase C8 as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.4 \sqrt{P}$	80 MHz to 800 MHz $d = 0.4 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.13	0.13	0.22
1	0.40	0.40	0.70
10	1.3	1.3	2.2
100	4.0	4.0	7.0

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Product / Produkt / Tuote / Product / Produkt / Προϊόν / Produkt

**bluephase C8 (G2)**

**GB** We herewith declare that the product listed above complies with the mentioned standards.  
Following the provisions of Directive(s):

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## Ivoclar Vivadent UK Limited
## Ivoclar Vivadent, Inc.